

Maintenance medium for human iPS cells

# ciKIC™ iPS medium

NEW!!

## *Feeder-free medium for human iPS cells*

The culture of human iPS cells (hiPSCs) requires frequent medium exchange, therefore the increase of workload and cost are seen as a serious problem. We developed new culture medium that can skip medium exchange on weekends and culture hiPSCs stably. ciKIC™ iPS medium supports your regenerative medicine research in the future.



### Advantages

- Low-protein & Albumin-free
- Weekend-free cell culture
- Single cell passaging
- High growth rate and stable pluripotency

### Product Information

Product name	Package	Storage	Product No.
ciKIC™ iPS medium	1kit (for 250 mL) Basal medium + Supplement Set	Basal medium : Store at 2 - 8°C Supplement Set : Store at -20°C	08371-13


- ★ This product was developed by technology transfer from Kyoto University, institute for integrated Cell-Material Sciences (iCeMS).
- ★ This product is for research use only. Do not use for human or animal diagnostic or therapeutic uses.

## Application


### Morphology

#### ■ Phase contrast image of iPS cells 7 days in culture

253G1




201B7



200  $\mu$ m

### Cell Growth

#### ■ Cell growth rate




**Normal hiPSCs morphology was observed in ciKIC™ iPS medium cultures, and ciKIC™ iPS medium supported stable high growth rate of hiPSCs.**


## Karyotype

#### ■ Karyotype analysis by Q-band

253G1



201B7



n=20 (253G1, 201B7)

**hiPSCs cultured in ciKIC™ iPS medium showed a normal karyotype.**


## Undifferentiation

### Flow cytometric analysis

\*Expression of undifferentiated marker was assessed by immunofluorescent stain and flow cytometer.

1<sup>st</sup> Ab ... 10 µg/mL anti TRA-1-60, anti SSEA-4

2<sup>nd</sup> Ab ... 2 µg/mL Alexa488 goat anti mouse IgG/IgM



hiPSCs cultured in ciKIC™ iPS medium expressed undifferentiated markers.


## Pluripotency

### Immunofluorescent staining images of iPSCs after differentiation

\*Embryoid bodies (EBs; using 253G1) were formed by floating culture. Then, EBs were transferred to gelatin-coated plate, continued culture and spontaneously differentiated.

1<sup>st</sup> Ab ... 5 µg/mL anti AFP, 1 µg/mL anti  $\alpha$ -SMA, 2 µg/mL anti  $\beta$ III-tubulin

2<sup>nd</sup> Ab ... 2 µg/mL Alexa488 goat anti mouse IgG/IgM



hiPSCs cultured in ciKIC™ iPS medium differentiated into three germ layers.

## Product Specification

### ▼ Basal medium



### ▼ Supplement Set



#### ■ Components

ciKIC™ iPS Basal medium : 1 bottle (250 mL)

ciKIC™ iPS medium Supplement Set

Supplement 1 : 1 tube (1 mL)

Supplement 2 : 1 tube (20 µL)

#### ■ Storage

ciKIC™ iPS Basal medium : Store at 2 - 8°C

ciKIC™ iPS medium Supplement Set : Store at -20°C

#### ■ Note

- Add Supplement 1 and 2 to Basal medium before use.
- Store at 2 - 8 °C in the attached aluminum package and use within 2 weeks after preparation of medium.

\* Contains Xeno-free components only.

\* Suitable for culturing on iMatrix-511 by pre-mix method or pre-coating method.

\* Enable single cell passaging by adding ROCK inhibitor (Y-27632).

\* Not need to add bFGF.

\* We recommend using ciKIC™ iPS medium for transfer from on feeder culture.

## Schedule of Medium Exchange

P : Passage, O : Medium exchange

	Thu.	Fri.	Sat.	Sun.	Mon.	Tue.	Wed.
Condition 1	P	O	—	—	—	O	O
Condition 2	P	O	—	—	O	O	O

⇒ Skip medium exchange on weekends (max 3 days in a row) .

Note: Select Condition 1 or Condition 2 depending on the cell line.

It confirmed that 253G1 and 201B7 can be cultured normally under Condition 1 and Condition 2.

## Reference

Yasuda, Shin-ya, et al., 2018, Chemically defined and growth-factor-free culture system for the expansion and derivation of human pluripotent stem cells. *Nature Biomedical Engineering* 2, 173-182

 **Kanto Chemical Co., Inc.**  
REAGENT DIVISION

2-1, Nihonbashi Muromachi 2-chome, Chuo-ku,  
Tokyo, 103-0022 JAPAN

E-mail : reag-info@gms.kanto.co.jp  
HP : <https://www.kanto.co.jp/english/>

May, 2019